

TEKNİK ÖZELLİKLER

İçerik : Örgü fiberglass ve su ile aktive olan poliüretan reçine

TECHNICAL PROPERTIES

Content : Knitted fiberglass and water-activated polyurethan resin.

EBAT / SIZES	EBAT (inç x yard) / SIZE (inches x yard)	KOLİ İÇİ (Adet) / BOX IN (pcs)	KOLİ EBATLARI (cm) / BOX SIZES (cm)
5 x 360	2" x 4	25	29 x 40 x 10
7,5 x 360	3" x 4	25	29 x 40 x 10
10 x 360	4" x 4	25	35 x 45 x 10
12,5 x 360	5" x 4	25	35 x 45 x 10